



Product Datasheet

Tie-2 (phospho-Tyr1108) antibody (orb770239)


Description	Rabbit polyclonal antibody to Tie-2 (phospho-Tyr1108)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human TIE2 around the phosphorylation site of Tyr1108. AA range:1074-1123
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent
Assay (Phosp...



Immunofluorescence
analysis of HepG2
cel...



Western blot
analysis of lysates
from NI...