



Product Datasheet

RGPD2/3/4/6/8 antibody (orb770344)

Description	Rabbit polyclonal antibody to RGPD2/3/4/6/8
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RGPD2/3/4/6/8. AA range:577-626
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110
Expiration Date	12 months from date of receipt.



Western blot analysis of the lysates from...



Western blot analysis of RGPD2/3/4/6/8 A...



Immunohistochemical analysis of paraffin...