


Product Datasheet

Vinculin (phospho-Tyr821) antibody (orb770365)


Description	Rabbit polyclonal antibody to Vinculin (phospho-Tyr821)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P
Immunogen	The antiserum was produced against synthesized peptide derived from human Vinculin around the phosphorylation site of Tyr821. AA range:786-835
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent Assay
(Phosp...



Immunofluorescence
analysis of HepG2
cel...



Immunohistochemistry
analysis of paraffi...