

Product Datasheet

Splicing factor 1 (phospho-Ser82) antibody (orb770411)

Description Rabbit polyclonal antibody to Splicing factor 1 (phospho-Ser82)

Species/Host Rabbit

Reactivity Human, Monkey, Mouse

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human SF1 around the phosphorylation site of Ser82. AA range:48-97

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 68


Expiration Date 12 months from date of receipt.



Western blot analysis of SH-SY5Y 293T us...



Immunofluorescence analysis of A549 cell...



Immunohistochemistry analysis of paraffi...