



Product Datasheet

CHSY2 antibody (orb770458)


Description	Rabbit polyclonal antibody to CHSY2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CHSY2. AA range:631-680
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of MCF7 cell...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from LO...