



Product Datasheet

BAIAP2L2 antibody (orb770486)


Description	Rabbit polyclonal antibody to BAIAP2L2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human BAIAP2L2. AA range:111-160
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from RA...



Western blot analysis of the lysates fro...



Immunohistochemical analysis of paraffin...