


Product Datasheet

Cubilin antibody (orb770490)

Description	Cubilin (CUBN) acts as a receptor for intrinsic factor-vitamin B12 complexes. The role of receptor i...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthesized peptide derived from the N-terminal region of human Cubilin
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	-20°C/1 year
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	Observed Band: 400kD
Dilution Range	Immunohistochemistry: 1/100 - 1/300
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...