

Product Datasheet

SMC1 (phospho-Ser966) antibody (orb770535)

Description Rabbit polyclonal antibody to SMC1 (phospho-Ser966)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human SMC1 around the phosphorylation site of Ser966. AA range:932-981

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 160


Expiration Date 12 months from date of receipt.



Western blot analysis of K562 using p-SM...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from HU...