



Product Datasheet

NIFK antibody (orb770629)


Description	Rabbit polyclonal antibody to NIFK
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human NIFK. AA range:200-249
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HUVEC cells...



Immunohistochemistry analysis of paraffin-embedded tissue sections...



Western blot analysis of lysates from K562 cells...