



Product Datasheet

LPAAT-Theta antibody (orb770655)


Description	Rabbit polyclonal antibody to LPAAT-Theta
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human AGPAT9. AA range:381-430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48
Expiration Date	12 months from date of receipt.



Western Blot analysis of HEPG2 using LPA...



Western blot analysis of SH-SY5Y 3T3 lys...



Immunofluorescence analysis of HepG2 cel...