



## Product Datasheet

### Mnk1 (phospho-Thr250) antibody (orb770699)

|                            |                                                                                                                                            |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Mnk1 (phospho-Thr250)                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                     |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                                                                                               |
| <b>Tested Applications</b> | ELISA, IHC-P                                                                                                                               |
| <b>Immunogen</b>           | The antiserum was produced against synthesized peptide derived from human MNK1 around the phosphorylation site of Thr250. AA range:216-265 |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4.                                                                                       |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                                    |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.                                                                                         |
| <b>Note</b>                | For research use only                                                                                                                      |
| <b>Isotype</b>             | IgG                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                            |



Enzyme-Linked  
Immunosorbent  
Assay (Phosp...



Immunohistochemical  
analysis of paraffin...