



Product Datasheet

PEA-15 (phospho-Ser116) antibody (orb770729)


Description	Rabbit polyclonal antibody to PEA-15 (phospho-Ser116)
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human PEA-15 around the phosphorylation site of Ser116. AA range:81-130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15
Expiration Date	12 months from date of receipt.



Western blot analysis of COLO 3T3 using ...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from CO...