


Product Datasheet

Actinin-alpha 1 antibody (orb770735)

Description	Rabbit polyclonal antibody to Actinin-alpha 1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Actinin- α 1 . at AA range: 1-80
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103
Expiration Date	12 months from date of receipt.



Western blot analysis
of 1) 3T3, 2) PC12...



Immunohistochemical
analysis of paraffin...