



Product Datasheet

RIP2 (phospho-Ser176) antibody (orb770748)


Description	Rabbit polyclonal antibody to RIP2 (phospho-Ser176)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RIPK2 around the phosphorylation site of Ser176. AA range:146-195
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Western Blot analysis of cell lysis treatment...



Western Blot analysis of 293 cells using...



Enzyme-Linked Immunosorbent Assay (Phospho...)