



Product Datasheet

RIP2 antibody (orb770749)


Description	Rabbit polyclonal antibody to RIP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RIPK2. AA range:146-195
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Western blot analysis of 3T3 lysis using...



Immunohistochemistry analysis of paraffi...



Western blot analysis of lysates from ra...