



Product Datasheet

Lfc antibody (orb770856)


Description	Rabbit polyclonal antibody to Lfc
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ARHGEF2. AA range:383-432
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	101
Expiration Date	12 months from date of receipt.



Western blot analysis of the lysates from...



Western blot analysis of lysates from RA...



Immunohistochemical analysis of paraffin...