



Product Datasheet

MSK1 (phospho-Thr581) antibody (orb770873)


Description	Rabbit polyclonal antibody to MSK1 (phospho-Thr581)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MSK1 around the phosphorylation site of Thr581. AA range:551-600
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90
Expiration Date	12 months from date of receipt.



Enzyme-Linked
Immunosorbent
Assay (Phosp...



Immunofluorescence
analysis of A549
cell...



Western blot
analysis of lysates
from RA...