



Product Datasheet

MPP1 antibody (orb770943)

Description	Rabbit polyclonal antibody to MPP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human KIF20B. AA range:1771-1820
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	210
Expiration Date	12 months from date of receipt.



Western blot analysis of KIF20B Antibody...



Immunohistochemical analysis of paraffin...