



Product Datasheet

CEP135 antibody (orb770957)


Description	Rabbit polyclonal antibody to CEP135
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CEP135. AA range:1081-1130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	134
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using CEP135 antibody.



Western blot analysis of mouse-lung and mouse-heart tissues using CEP135 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CEP135 antibody.