

Product Datasheet

Histone H3 (Mono Methyl Lys27) antibody (orb771040)

Description	Rabbit polyclonal antibody to Histone H3 (Mono Methyl Lys27)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide of Histone H3 (Mono Methyl Lys27)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15-17
Expiration Date	12 months from date of receipt.