

Product Datasheet

HSC 70 antibody (orb771061)

Description	Rabbit polyclonal antibody to HSC 70
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthesized peptide derived from the Internal region of human HSC 70.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70-72
Expiration Date	12 months from date of receipt.