

## Product Datasheet

### PPAR-gamma antibody (orb771132)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to PPAR-gamma             |
| <b>Species/Host</b>        | Rabbit                                               |
| <b>Reactivity</b>          | Human, Mouse, Rat                                    |
| <b>Conjugation</b>         | Unconjugated                                         |
| <b>Tested Applications</b> | IHC-P, WB                                            |
| <b>Immunogen</b>           | Recombinant Protein of PPAR-γ                        |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze-thaw cycles.   |
| <b>Note</b>                | For research use only                                |
| <b>Isotype</b>             | IgG                                                  |
| <b>Clonality</b>           | Polyclonal                                           |
| <b>MW</b>                  | 53,57                                                |
| <b>Expiration Date</b>     | 12 months from date of receipt.                      |