

Product Datasheet

Actinin-alpha 1 antibody (orb771134)

Description	Rabbit polyclonal antibody to Actinin-alpha 1
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein of Actinin- α 1
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103
Expiration Date	12 months from date of receipt.