



Product Datasheet

CD26 antibody (orb771262)


Description	Rabbit polyclonal antibody to CD26
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human DPP4. AA range:350-400
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85
Expiration Date	12 months from date of receipt.



Western blot analysis of 293T 3T3 lysis ...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...