



Product Datasheet

CD249 antibody (orb771288)


Description	Rabbit polyclonal antibody to CD249
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human ENPEP. AA range:671-720
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110
Expiration Date	12 months from date of receipt.



Western Blot analysis of MCF7 cells using CD249 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CD249 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CD249 antibody.