



Product Datasheet

C7 antibody (orb771297)


Description	Rabbit polyclonal antibody to C7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human C7. AA range:100-150
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93
Expiration Date	12 months from date of receipt.



Western blot analysis of Hela lysis using C7 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using C7 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using C7 antibody.