

Product Datasheet

Acetyl-EPAS-1 (K385) antibody (orb771320)

Description	Rabbit polyclonal antibody to Acetyl-EPAS-1 (K385)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized acetyl-peptide of Acetyl-EPAS-1 (K385)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.