



Product Datasheet

Perforin 1 antibody (orb771385)


Description	Rabbit polyclonal antibody to Perforin 1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human PRF1. AA range:451-500
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	61
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using Perforin 1 antibody (orb771385).



Immunohistochemical analysis of paraffin-embedded tissue sections using Perforin 1 antibody (orb771385).



Immunohistochemical analysis of paraffin-embedded tissue sections using Perforin 1 antibody (orb771385).