



Product Datasheet

TyrRS antibody (orb771391)


Description	Rabbit polyclonal antibody to TyrRS
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human YARS. AA range:451-500
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using TyrRS antibody (orb771391).



Immunohistochemical analysis of paraffin-embedded tissue sections using TyrRS antibody (orb771391).



Immunohistochemical analysis of paraffin-embedded tissue sections using TyrRS antibody (orb771391).