


Product Datasheet


GRP78 BiP antibody (orb771482)

Description	Rabbit polyclonal antibody to GRP78 BiP
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide from human protein at AA range: 505-570
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	78
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1,

...



Western blot analysis of 3T3 Cell Lysate...



Immunohistochemical analysis of paraffin...