



Product Datasheet

Calcitonin antibody (orb771484)


Description	Rabbit polyclonal antibody to Calcitonin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthetic peptide from human protein at AA range: 80-141
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat U937 Raj...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...