



Product Datasheet

JAK2/3 (Phospho-Tyr966/939) antibody (orb771687)


Description	Rabbit polyclonal antibody to JAK2/3 (Phospho-Tyr966/939)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized Phospho peptide derived from human JAK2/3.at AA range: T966
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western Blot analysis of 1, mouse-liver ...



Western Blot analysis of hela cells usin...



Immunohistochemical analysis of paraffin...