


Product Datasheet

ACSS1 (Acetyl-K642) antibody (orb771693)

Description	Rabbit polyclonal antibody to ACSS1 (Acetyl-K642)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized Acetyl peptide derived from human ACSS1. at AA range: K642
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western
Blot analysis
of 1, 293T
2, hela...