



Product Datasheet

Smad3 (Phospho-Ser213) antibody (orb771752)

Description	Rabbit polyclonal antibody to Smad3 (Phospho-Ser213)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho derived from human Smad3 (Phospho-Ser213)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western blot analysis of various lysate,...



Immunohistochemical analysis of paraffin...