



Product Datasheet

TIE2 antibody (orb771767)


Description	Rabbit polyclonal antibody to TIE2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from human TIE2 Polyclonal
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.



Western blot analysis
of CACO2 lysate, a...



Immunohistochemical
analysis of paraffin...



Immunohistochemical
analysis of paraffin...