


Product Datasheet

PNISR antibody (orb773186)

Description	Rabbit polyclonal antibody to PNISR
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from part region of human protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	88
Expiration Date	12 months from date of receipt.



Western Blot analysis of Hela lysis, using PNISR antibody.