

Product Datasheet

Mouse ACE ELISA Kit (orb774866)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.32-20 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cerebrospinal fluid, cell culture supernates and other biological fluids
Uniprot ID	P09470
Sensitivity	0.142 ng/mL
Expiration Date	Please enquire.

