

Product Datasheet

Human CLU ELISA Kit (orb774944)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	12.5-800 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, urine, saliva, tissue homogenates, cerebrospinal fluid, cell culture supernates and other biological fluids
Uniprot ID	075153
Sensitivity	5.1 ng/mL
Expiration Date	Please enquire.

