

## Product Datasheet

### Mouse NGAL ELISA Kit (orb774950)

|                            |                                                                                                             |
|----------------------------|-------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...     |
| <b>Reactivity</b>          | Mouse                                                                                                       |
| <b>Range</b>               | 0.16-10 ng/mL                                                                                               |
| <b>Tested Applications</b> | ELISA                                                                                                       |
| <b>Note</b>                | For research use only                                                                                       |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, urine, cell culture supernates and other biological fluids |
| <b>Uniprot ID</b>          | <b>P11672</b>                                                                                               |
| <b>Sensitivity</b>         | 0.057 ng/mL                                                                                                 |
| <b>Expiration Date</b>     | Please enquire.                                                                                             |

