

Product Datasheet

Human APOA1 ELISA Kit (orb774957)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	6.25-400 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P02647
Sensitivity	2.73 ng/mL
Expiration Date	Please enquire.

