

Product Datasheet

Mouse MIP1b ELISA Kit (orb774983)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	12.5-800 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P14097
Sensitivity	5.3 pg/mL
Expiration Date	Please enquire.

