

Product Datasheet

Human LEP ELISA Kit (orb775025)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, saliva, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P41159
Sensitivity	0.058 ng/mL
Expiration Date	Please enquire.

