

Product Datasheet

Rabbit TNFa ELISA Kit (orb775029)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Rabbit
Range	3.13-200 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P04924
Sensitivity	1.17 pg/mL
Expiration Date	Please enquire.

