

Product Datasheet

Mouse Kim1 ELISA Kit (orb775030)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	39.07-2500 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluidstissue homogenates
Uniprot ID	Q5QNS5
Sensitivity	15 pg/mL
Expiration Date	Please enquire.

