

Product Datasheet

Mouse TGFb1 ELISA Kit (orb775051)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	15.63-1000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids
Uniprot ID	P04202
Sensitivity	6.5 pg/mL
Expiration Date	Please enquire.

