

Product Datasheet

Human VEGFC ELISA Kit (orb775059)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	78.13-5000 pg/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P49767
Sensitivity	33 pg/mL
Expiration Date	Please enquire.

