

Product Datasheet

Human VCAM1 ELISA Kit (orb775121)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Human
Range	0.79-50 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P19320
Sensitivity	0.29 ng/mL
Expiration Date	Please enquire.

