

Product Datasheet

Mouse TIMP1 ELISA Kit (orb775125)

Description	The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided...
Reactivity	Mouse
Range	0.16-10 ng/mL
Tested Applications	ELISA
Note	For research use only
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids
Uniprot ID	P12032
Sensitivity	0.059 ng/mL
Expiration Date	Please enquire.

