

## Product Datasheet

### Guinea pig IL1b ELISA Kit (orb775146)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Guinea pig                                                                                              |
| <b>Range</b>               | 7.82-500 pg/mL                                                                                          |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids    |
| <b>Uniprot ID</b>          | <a href="#">Q9WVG1</a>                                                                                  |
| <b>Sensitivity</b>         | 2.8 pg/mL                                                                                               |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |

