

## Product Datasheet

### Human IL20 ELISA Kit (orb775156)

|                            |                                                                                                         |
|----------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | The test principle applied in this kit is Sandwich enzyme immunoassay. The microtiter plate provided... |
| <b>Reactivity</b>          | Human                                                                                                   |
| <b>Range</b>               | 15.63-1000 pg/mL                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                   |
| <b>Note</b>                | For research use only                                                                                   |
| <b>Sample Types</b>        | Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Uniprot ID</b>          | <b>Q9NYY1</b>                                                                                           |
| <b>Sensitivity</b>         | 6 pg/mL                                                                                                 |
| <b>Expiration Date</b>     | Please enquire.                                                                                         |

